eng TwoFer™ Huber-Tipped Needle

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

ita Ago con punta Huber TwoFer


spa Aguja de punta Huber TwoFer

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swe TwoFer Huber-nål

<div> <div> <div></div> <div></div> </div> <div> <div>CE </div> <div>2797 </div> </div> </div>	<p>MedSafety Solutions 7012 S. Revere Pkwy, Suite 120 Centennial, CO 80112 USA www.medsafetysolutions.com</p>	<div> <div> <div></div> <div></div> </div> <div> <div>EU </div> <div>EMERGO EUROPE </div> </div> </div> <p>Prinsessegrocht 20 2514 AP The Hague The Netherlands</p>
© Copyright 2021 MedSafety Solutions		53001247 Rev. F (07/2021)
<p>1.</p> 		

2.		
----	---	--

3.		
----	--	--

ita Ago con punta Huber TwoFer

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swe TwoFer Huber-Tipped Needle

Available in both long (green) and short (purple) versions, these 16-gauge thin-walled needles function as both fluid transfer devices and as venting needles. Hubs are colorcoded to assist with visual identification.

Specifications

- Transfers liquids into closed vials
- Transfers high-volume liquids using pumps and repeating syringes
- Sterile and non-pyrogenic in unopened, undamaged cap and sheath (sterile barrier)
- Length:
 - Long – approximately 1” cannula plus 0.5” vent/hub
 - Short – approximately 0.5” cannula plus 0.5” vent/hub
- Huber Point
- Not made with natural rubber latex
- Contains no di(2-ethylhexyl) phthalate additives
- Made in the USA

Materials

- Cannula: 304 Stainless Steel
- Hub/Vent/Cap/Sheath: Polypropylene
- Adhesive: Acrylated Urethane
- Lubricant: Medical Grade Silicone

Indications for Use

Designed as a needle for withdrawing medication from a vial using the metal cannula and for adding to a vial by using the plastic vent hub to equalize pressure.

Contraindications: For Pharmacy Use Only, not for use with Blood Administration.

Precautions

- Venting opens the container to the environment.
- Inspect all compounded products for particulate, per government/clinical recommendations.
- Verify compatibility between drugs to be transferred and materials used in the manufacturing of this product with the drug manufacturer’s package insert and/or material safety documentation.
- When preparing compounded sterile products, aseptic technique must be maintained throughout the production process.
- Avoid long-term exposure to direct light.
- Do not use product if protective cap has fallen off or is damaged.
- Do not reuse. Product reuse can be unsafe because of the risk of infection, cross contamination, endotoxin reaction, patient injury due to material failure and/or chemical reaction due to cleaning or sterilization agents.
- Dispose of properly according to applicable sharps disposal rules and regulations.

Instructions for Use

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and appropriate regulatory authorities. Within the EU, report to the competent authority of the Member State in which the user and/or patient is established.
1. Twist off cap and attach luer connection to syringe or tube set. Remove sheath.

- For vented applications:
 - Insert both the stainless steel needle cannula and plastic venting hub through the septum.
 - Note - If holding vial in an inverted position, do not attempt to expel air or fluid from the syringe with vent inserted.
- For non-vented applications:
 - Insert only the stainless steel needle cannula portion of the TwoFer Needle through the septum.

SYMBOL GLOSSARY	
Symbol	Title of Symbol / Definition / Reference
	Do not re-use / Indicates device is intended for one use or for use on a single patient during a single procedure / Ref. 5.4.2 ISO 15223-1:2016
	Consult instructions for use / Indicates the need for the user to consult the instructions for use / Ref. 5.4.3 ISO 15223-1:2016
	CAUTION! / Indicates important cautionary information. / Ref. 5.4.4 ISO 15223-1:2016
	Do not use if package is damaged / Indicates device should not be used if package is damaged or opened / Ref. 5.2.8 ISO 15223-1:2016
	Manufacturer / Indicates the medical device manufacturer name & address / Ref. 5.1.1 ISO 15223-1:2016
 LOT	Batch code / Indicates the manufacturer’s batch code / Ref. 5.1.5 ISO 15223-1:2016
 REF	Catalogue Number / Indicates the manufacturer’s catalogue number / Ref. 5.1.6 ISO 15223-1:2016
	Use-by-date / Indicates the date after which the device is not to be used / Ref. 5.1.4 ISO 15223-1:2016
 STERILE IN	Sterilized using irradiation / Indicates device has been sterilized using irradiation / Ref. 5.2.4 ISO 15223-1:2016
 Rx Only	Prescription Device / CAUTION! Federal (US) law restricts this device to sale, distribution, and use by or on order of a physician / Ref. 21 CFR § 801.109
 MD	Medical Device / Indicates product is a Medical Device / Ref. MDR 2017/745 Annex MD, Chapter III, 23.2 (g)

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